

Technical data sheet
Product groups 08 Cannula retaining strap



Hook- & loop tape

Hook tape	ÖKO-TEX Standard 100 Class II 100% Polyamide PA
-----------	--

Loop tape	ÖKO-TEX Standard 100 Class II 100% Polyamide PA
-----------	--

Soft foam

Foam RG26 lined in both sides with jersey and velor. Manufactured according to ÖKO-TEX Standard 100 Class II. Washable up to 40°C. Minor shrinkage, similar to cotton, possible.

Rubber band

Manufactured according to ÖKO-TEX Standard 100 Class II