

Technical data sheet

Product group 03 universal tube holder



Hook- & loop tape

Hook tape	ÖKO-TEX Standard 100 Class II 100% Polyamide PA
Loop tape	ÖKO-TEX Standard 100 Class II 100% Polyamide PA
Chemical resistance	resistant to - diluted acids and alkaline chemicals - Trichlorethylene - solvent at chem. cleaning
Cleaning instructions	- was hup to 60°C - chemical cleaning - dry-spin - sterilize in autoclave

Autoclaving at 100°C to 134°C will attack the hooks, which will affect adhesion after several processes. In any case, keep the hook and loop tape closed.

CLIP with spring / HOOK

Consisting of stainless steel.

Rivets

Consisting of nickel-plated brass, rustproof

Manufacturer

GERWING Medizinprodukte
Gerd Riester e.K.
Hochmeisterstr. 55
D-72417 Jungingen